

In Unison™



In Unison™

Using SoftView and OnGuard together offers exceptional visualization of chest X-rays, without interrupting your workflow. SoftView analyzes the chest X-ray to produce a bone suppressed soft tissue image for improved clarity while OnGuard identifies solitary pulmonary nodules that may be lung cancer.¹ The technologies are automatically presented “in unison” as a SoftView image with OnGuard CAD markings, without interrupting workflow. Displaying this integrated image is as easy as the click of a mouse.




SoftView Enhanced Chest Imaging is an innovative image enhancement technology designed to suppress ribs and clavicles on chest X-rays, creating a soft tissue image for every chest X-ray, even those taken on portables. SoftView requires no additional radiation dose or patient procedure and aids radiologists in the detection of more lung nodules.



OnGuard Computer Aided Detection (CAD) identifies solitary pulmonary nodules that may be early stage lung cancer. The OnGuard algorithm automatically processes the chest X-ray and places scaled markers over regions of interest to improve the detection of 9-30mm nodules.

Detect more with SoftView and OnGuard.

Using SoftView Enhanced Chest Imaging and OnGuard Chest X-ray CAD is as simple as a click of the mouse. Improving the clarity of soft tissue and detecting lung cancer in its earliest stages has never been easier.



Standard Digital X-ray



SoftView Enhanced Image



SoftView/OnGuard Combined




1 <http://lungcancer.about.com/od/typesoflungcancer/a/spn.htm>

➤ Truth Confirmed – Lung Nodule

Riverain Medical

Riverain Medical believes in doing more with less. Our technologies are designed to enhance the expert skills of radiologists using existing equipment to improve patient outcomes without additional radiation dose or procedures. We are driven to save lives by improving the detection and diagnosis of diseases and other clinical conditions. As an industry leader and technology innovator, our practical solutions provide physicians with tools to assist in image interpretation



To learn more:

Contact us at (800) 990-3387
or visit
www.riverainmedical.com.



Riverain Medical
3020 South Tech Boulevard
Miamisburg, Ohio 45342-4860
P: +1 800.990.3387 +1 937.425.6811
F: +1 937.425.6493