


Equipment Layout



RECOMMENDED CEILING HEIGHT 9'-0"

07159-A1

EQUIPMENT LEGEND

CFCI.....CONTRACTOR FURNISHED AND CONTRACTOR INSTALLED
 OFCI.....OWNER FURNISHED AND OWNER OR CONTRACTOR INSTALLED
 OFOI.....OWNER FURNISHED AND OWNER OR VENDOR INSTALLED

	ITEM	DESCRIPTION	SIZE	WT/LBS	BTU/HR	DETAIL SHEET
OFOI	①	QUANTUM ODYSSEY QG-65 CONSOLE ON SHELF BELOW COUNTER	18" W 12" D 3" H	5	N/A	
OFOI	②	QUANTUM ODYSSEY QG-65 GENERATOR CABINET	24" W 20" D 34" H	380	4000	
OFOI	③	QUANTUM QT-750 RADIOGRAPHIC TABLE	85" W 36" D 33" H	500	230	
OFOI	④	QUANTUM QS-550 FLOOR MOUNTED TUBESTAND	108" L 48" D 92" H	600	873	
OFOI	⑤	QUANTUM QW-420 VERTI-Q WALL STAND	26" W 23" D 84" H	225	N/A	
OFOI	⑥	CANON DR SYSTEM CONTROL PANEL/KEYBOARD	18" W 8" D 5" H	6	100	
OFOI	⑦	CANON DR SYSTEM C3 CONTROL CABINET (UNDER COUNTER)	12" W 21" D 23.4 H	77	854	
OFOI	⑧	CANON CXDI SENSOR POWER SUPPLY UNIT	4.3" W 15.7" D 6.3" H	10	100	
OFOI	⑨	CANON DR SYSTEM MULTIBOX / HUB (UNDER COUNTER)	3.9" W 13" D 8.5" H	9	100	
CFCI	⑩	LEAD GLASS VIEW WINDOW	24" W 36" H	N/A	N/A	
CFCI	⑪	COUNTER TOP (30-36" AFF)	24" D	N/A	N/A	

CUSTOMER:

Quantum DS-3 Radiographic System Installation

THESE DRAWINGS ARE SUBMITTED TO SUGGEST LOCATION OF ALL DIAGNOSTIC IMAGING EQUIPMENT AND ASSOCIATED APPARATUS. ELECTRICAL WIRING DETAILS AND ROOM ARRANGEMENTS. IN PREPARING THIS DRAWING, EVERY EFFORT HAS BEEN MADE TO CONFORM TO THE ACTUAL EQUIPMENT EXPECTED TO BE INSTALLED. HOWEVER, THESE DRAWINGS ARE NOT TO BE USED FOR ACTUAL CONSTRUCTION PURPOSES, AND THE PREPARER CANNOT ACCEPT ANY RESPONSIBILITY FOR DAMAGES RESULTING THEREFROM. CONSTRUCTION DRAWINGS TO BE PREPARED BY PURCHASER'S ARCHITECT OR CONTRACTOR.

CMS IMAGING, INC.
 4050 Azalea Drive
 North Charleston, SC 29405 (843) 763-1016

DRAWN BY:
R. Rider
 DATE:
05/17/07
 DWG. NO.
07159
 SHEET:
A1