


Control

NOTE: WINDOW BEHIND WALL STAND
TO BE REMOVED AND WALL
BLOCKED UP

Equipment Layout



EXISTING CEILING HEIGHT 9'-3"
RECOMMENDED CEILING HEIGHT 9'-0"

08261A-A1

EQUIPMENT LEGEND

CFCI.....CONTRACTOR FURNISHED AND CONTRACTOR INSTALLED
OFCI.....OWNER FURNISHED AND OWNER OR CONTRACTOR INSTALLED
OFOI.....OWNER FURNISHED AND OWNER OR VENDOR INSTALLED

	ITEM	DESCRIPTION	SIZE	WT/LBS	BTU/HR	DETAIL SHEET
OFOI	①	QUANTUM ODYSSEY SERIES CONSOLE	18" W 12" D 3" H	5	N/A	
OFOI	②	QUANTUM ODYSSEY GENERATOR CABINET	24" W 20" D 34" H	380	4000	
OFOI	③	QUANTUM QT-750 RADIOGRAPHIC TABLE	79" W 36" D 33" H	500	230	
OFOI	④	QUANTUM RS-580 FIXED CEILING RAILS	135" L	270	N/A	
OFOI	⑤	QUANTUM RS-580 MOVING TUBE CARRIAGE	119" W 24" D 45" H	796	873	
OFOI	⑥	QUANTUM QW-420 VERTICAL WALL STAND	24" W 14" D 86" H	130	N/A	
CFCI	⑦	EXISTING LEAD GLASS VIEW WINDOW	N/A	N/A	N/A	
CFCI	⑧	COUNTER TOP (36-40" AFF)	18" D	N/A	N/A	

CUSTOMER:

Equipment Layout / Legend

THESE DRAWINGS ARE SUBMITTED TO SUGGEST LOCATION OF ALL DIAGNOSTIC IMAGING EQUIPMENT AND ASSOCIATED APPARATUS, ELECTRICAL WIRING DETAILS, AND ROOM ARRANGEMENTS. IN PREPARING THIS DRAWING, EVERY EFFORT HAS BEEN MADE TO CONFORM DETAILS TO THE ACTUAL EQUIPMENT EXPECTED TO BE INSTALLED. HOWEVER, THESE DRAWINGS ARE NOT TO BE USED FOR ACTUAL CONSTRUCTION PURPOSES, AND THE PREPARER CANNOT ACCEPT ANY RESPONSIBILITY FOR DAMAGES RESULTING THEREFROM. CONSTRUCTION DRAWINGS TO BE PREPARED BY PURCHASER'S ARCHITECT OR CONTRACTOR.



DRAWN BY:
R. Rider
DATE:
08/18/08
DWG. NO.
08261A
SHEET:
A1