




interWORKS *solutions for* Orthopedics

Is your orthopedic PACS living up to your expectations?

Orthopedic surgeons face an ever-expanding set of challenges each day, including the pressure to reduce costs and deliver faster, higher-quality patient care. To meet these challenges, Avreo designed interWORKS as a WAN optimized, web-based system, that streamlines your orthopedic imaging workflow in the office, operating room, and wherever your practice takes you. Its zero-impact

installation design does not require third-party utilities, which allows it to go where many viewers can't. Its simple, yet powerful viewer automatically displays images according to each orthopedic surgeon's preferences. In addition, mobile reference tools provide powerful visual illustrations when discussing treatment plans with patients.

EHR for Radiology because RIS/PACS is not enough



What Do I Need to Know?

Our organization uses a radiology group to provide interpretation for specific studies. Can you send these studies to their PACS?

Yes. A technologist can send these studies as needed to your radiology reading group. In addition, your radiology group can read remotely from interWORKS if needed.

Can interWORKS receive orders from my Practice Management/EHR?

Yes. To maximize workflow efficiency, interWORKS can receive orders from your practice management/EHR system via HL7. Our visual DayPlanner then displays these orders to help technologists manage their schedule within a single application.

Can interWORKS provide radiology images within my Practice Management/EHR?

Yes. interWORKS supports this functionality to assist your organization in meeting the Image Results Meaningful Use measure.

Do you support orthopedic templating tools?

Yes. Avreo partners with OrthoView to provide Orthopedic templating. OrthoView enables surgeons to plan and template orthopaedic procedures with digital images from interWORKS and reduce or eliminate their need for printed film. Templated images and the OrthoView planning report are available in interWORKS for future reference.

interWORKS solutions for

Orthopedics

Recommended Configuration

- interWORKS
- interWORKS VM hosted in your organization's data center. This option requires Avreo approved VMWare or HyperV virtual guest on your organization's hardware.
- SLS or ASP licensing.

Why interWORKS?


- Easily accepts radiology orders from your practice management/EHR system to minimize impact on your scheduling workflow.
- Includes quality assurance tools for technologists to review image quality and order accuracy before submitting a study to the orthopedic surgeon.
- Empower orthopedic surgeons to read virtually anywhere. interWORKS offers WAN optimized, web-based image and reporting access to ensure all healthcare professionals can do their job regardless of location.
- Status driven workflow configurable to the needs of orthopedic practices.
- Quickly document case notes via digital dictation or quick-pick templates.
- Storage agnostic and vendor neutral for archiving. Avreo has experience connecting with multiple VNAs.

Proven Quality Over Time

Your healthcare organization makes a commitment to your patients to provide high quality medical care in a safe and welcoming environment. We would like the opportunity to demonstrate this same level of excellence and commitment to your healthcare organization.

Since 2008, Avreo has been named KLAS Category Leader in PACS/Radiology 6 times.

Easily Access Radiology Images and Reports Throughout Your Enterprise



Let's Start a Conversation!

We want to discuss the current challenges that you are facing so we can fully understand the needs of your organization and offer solutions to alleviate these issues.

To Learn More

 www.avreo.com

 1-800-844-9926

 Live demonstrations

