

See More. See it Fast. See it First.

Konica Minolta's AeroDR[®] HD by the Numbers.*



AI READY

4X MORE DATA IN EVERY IMAGE

HD 3488 x 4256: 14.8 M



10% HIGHER X-RAY SENSITIVITY

XE = 65% (@ 0 cycle/mm)
HD = 72% (@ 0 cycle/mm)



7 YEAR POWER CELL WARRANTY

No batteries to change



UP TO 8 HOURS OF UPTIME ON A 30 MINUTE CHARGE



EXPOSURE LENGTHS UP TO 4 SECONDS

with Automatic Exposure Detection
(10 seconds with generator connection)



REALISM[™] ADVANCED HD IMAGE PROCESSING



Built to last



BEND RESISTANCE
Up to 286 lbs



20% HIGHER DROP RESISTANCE
Over 1.2 Meters



WATERPROOF
IPX6 Certified



STRONGER
20% Higher Point Load
33% Higher Distributed Load



EASIER TO HOLD
60% Less Grip Strength Needed

* All specifications may vary by size, installation and usage conditions. Please refer to current product manual for latest specifications.