

DEL MEDICAL

MT500

Mobile Table

The Del Medical MT500 mobile radiographic table is ideally suited for U-Arm systems located in imaging centers, orthopedic clinics and private practices. It can also be paired with an overhead tube crane or floor mounted tube stand to offer an economical table option. It features a heavy-duty steel frame for maximum stability and a stain-resistant fiber-resin tabletop. The large-diameter casters provide smooth movement, swivel 360 degrees, and have easy, yet reliable, foot locks.

Product Features

- All steel structure and welded base provides maximum strength and reliability
- 4 durable swivel casters with foot locks
- Smooth, stain resistant fiber resin top
- 500 lb. (227kg) patient weight limit



MT500

Mobile Table



Technical Specifications

- Table top length: 80.0" (203.2cm)
- Table top width: 26.0" (66.1cm)
- Table top height: 27.9" (70.9cm)
- Weight without cabinet: 164 lbs. (74kg)
- Weight with cabinet: 201 lbs. (91.2kg)






Table Options

- Grid cabinet with tray
 - Focal spot travel length: 43.0" (109.22cm)
- Table side patient handgrips
- Lateral cassette holder



* All specifications subject to change without notice



02.17 #9000-MT500

DEL MEDICAL

Toll Free: (800) 800-6006
Fax: (847) 288-7011
Email: info@delmedical.com
www.delmedical.com
241 Covington Drive, Bloomingdale, IL 60108

UMG

Toll Free: (800) 261-9808
Fax: (914) 835-6111
Email: info@umgxr.com
www.umgxr.com
28 Calvert Street, Harrison, NY 10528