

# RT 100

## FOUR-WAY FLOAT TOP MEDICAL X-RAY TABLE

This sturdy table is ideally suited for a variety of general radiographic needs including small hospitals, clinics, ambulatory care centers and private offices. The floating table top glides freely on a precision roller bearing system and can support patients up to 700 lbs. (318 kg).



*Foot treadle electric lock release runs  
the entire length of pedestal.*

### Product Features

- 17" x 17" image receptor with heavy-duty manual tray
- Heavy-duty electric table locks with foot treadle release
- All steel structure and welded base provides maximum strength and reliability
- Smooth, stain resistant fiber resin top
- Pedestal base can house generator electronics to save space

