
VS100


VERTICAL WALL STAND

A counterbalanced wall stand with vertical electric locks that allow for easy positioning. The extensive vertical travel allows for a complete range of exams from skull to weight bearing knee and ankle exams.

Product Features

- Dual counterweight cables for added safety
- Vertical travel range of 58.38" (148.28cm)
- Available with 17" x 17" (43cm x 43cm) grid cabinet or bucky and a wide slection of grid options







VS100

Vertical Wall Stand

Technical Specifications

Wall stand height: 86.75" (220.4cm)
Receptor height: 23.38" (59.4cm)
Vertical travel range: 58.38" (148.28cm)
Wall stand width (with receptor): 23.25" (59cm)
Wall stand depth (with receptor): 13" (33cm)

Maximum focal spot height: 72.38" (183.84cm)

Minimum focal spot height: 14" (35.6cm)

Front panel aluminum

equivalency: 0.4mm

Weight: 258lb (117kg)

Power requirements: 24VDC




Certifications

ETL Listed

Product Options

- 17" x 17" (43cm x 43cm) grid cabinet or bucky
- Deluxe, heavy-duty manual cassette tray
- Size-sensing cassette tray with PBL interface
- Ion chamber (AEC)
- · Variety of grid options
- 24VDC external power supply 110/220V

Optional cassette tray accepts 17" x 17" (43cm x 43 cm) or smaller cassettes or digital panels.



DEL MEDICAL

Toll Free: (800) 800-6006 Fax: (847) 288-7011 Email: info@delmedical.com www.delmedical.com Del Medical is a UMG Company. Del Medical is an ISO 13485 & ISO9001 Certified Facility. Specifications are subject to change without notice. UMG

Toll Free: (800) 261-9808 Fax: (914) 835-6111 Email: info@umgxray.com www.umgxray.com 28 Calvert Street, Harrison, NY 10528