


CL24s 24 Inch LED DISPLAY



i Product dimension



i Display Performance

Model		CL24s
LCD Panel	Description	24.1" TFT LCD
	Resolution	1920 x 1200 pixel
	Number of colors	16.7 M colors
Visual Performance	Brightness	300 cd/m ² (Typical)
	Contrast Ratio	1000 : 1 (Typical)
	Viewing Angle	R/L 178°, U/D 178° where CR < 10
Input Signal		2x DVI-D
		1 x VGA
		1 x C-VIDEO
		1 x S-VIDEO
External Control		RS-232C (Mini USB JACK)
Power Supply		AC/DC Adaptor (AC 100-240V~, DC 24V 3.75A)
Power Consumption		60W (Max)
Control Key		Power, Menu, PIP, Up, Down, Plus, Minus, Input
Active Area		518.4mm(W) X 324.0mm(H)
Unit Dimension		565.9mm(W) x 460.7mm(H) x 224.5mm(D)
Weight		9.42 kg / 20.77 lbs

i Environmental Consideration

Operating Conditions	Temperature	0° ~ 40°
	Humidity	5% ~ 85%, non-condensing
Transportation Conditions	Temperature	-20° ~ 60°
	Humidity	5% ~ 95%, non-condensing
Storage Conditions	Temperature	-20° ~ 60°
	Humidity	30% ~ 80%, non-condensing
Regulatory Approvals	EMC	KC : KN22, KN24 FCC : Part 15 Class B CE-EMC : EN55011, EN60601-1-2, EN61000-3-2, EN61000-3-3
	ITE Safety	CB : IEC(EN)60950-1 2 _{nd} Edition
	Medical Safety	UL/cUL : UL60601-1, CAN/CSA-C22.2 No.601.1-M90 CB : IEC(EN)60601-1 3 _{rd} Edition
	CCC	GB4943.1, GB9254, GB17625.1
	CEL	GB21520



WIDE Corporation (Gomae-Dong) 12 Wongomae-Ro, Giheung-Gu, Yongin-Si, Gyeonggi-Do 446-901, Korea Tel: +82-31-218-1600 Fax: +82-31-274-7400 info@widecorp.com **WIDE USA Corporation** 2210 E. Winston Road, Anaheim, CA 92806, USA Tel: +1-714-300-0540 infousa@widecorp.com **WIDE EUROPE B.V.** Panamalaan 3L, 1019 AS Amsterdam, The Netherlands Tel: +31-20-311-9797 Fax: +31-20-311-9790 infoeu@widecorp.com **WIDE JAPAN** Shinjuku Suzuki Bldg A 1-6-8 Shinjuku, Shinjuku-Ku, Tokyo 160-0022, Japan Tel: +81-3-6457-8371 Fax: +81-3-6457-8372 infojapan@widecorp.com **WIDE ASIA** 1350/88-89, Thairong tower, Pattanakarn road, Suanluang, Bangkok 10250, Thailand Tel: +66 2 717 1843 Fax: +66 2 717 1844 infoasia@widecorp.com **Foreseeson Custom Displays, Inc.** 2210 E. Winston Road Anaheim, CA 92806 USA Tel: +1 714 300 0540 Fax: +1-714-300 0546 **Foreseeson GmbH** Industriestrasse 38a 63150 Heusenstamm, Germany Tel: +49 6104 64398 0 Fax: +49 6104 64398 11 sales@fsnmed.eu **Foreseeson UK Ltd.** Unit 2 Kingsmill Business Park Chapel Mill Road Kingston upon Thames Surrey KT1 3GZ Tel: +44 208 546 1047 Fax: +44 208 546 3931 sales@fsnmed.eu | Created in Korea. 2014. 08. ver 1.0