




MX50N 21.3" Mono Display

- 2048 x 2560
- 1,200cd/m²
- Embedded IQ Sensor[®]
- DisplayPort & DVI



i Display Performance

| | |
|---------------------|--------------------------------------------------------------|
| Screen Size | 21.3" |
| Viewable Image Size | 422.4mm x 337.92mm (16.6"x13.3") |
| Display Resolution | 2048(H) x 2560(V) |
| Brightness | 1,200 cd/m ² |
| Contrast Ratio | 1000:1 |
| Viewing Angle | 89° / 89° / 89° / 89° (U/D/L/R) |
| Input Signal | DVI-D, Display Port |
| Control Key | Tact button |
| Power Supply | AC Input : 100-240Volt±10% 60Hz/50Hz±3Hz DC Output : +24V |
| Power Consumption | 50W |
| Power Save Mode | 2W |

| | |
|------------------|--------------------------------------|
| Backlight System | LED |
| Built-in Sensors | IQ Sensor 3 rd , SBC, DAC |
| Dimension | 390.3mm(W)x520.1mm(H)x248.8mm(D) |
| Weight | 10.7Kg(23.59lb) with Stand |
| Options | Protective Glass |

i Environmental Consideration

| | |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| Temperature | Operating : 0°C ~ 40°C |
| | Storage : -20°C ~ 60°C |
| Humidity | Operating : 10% ~ 80% |
| Certifications and Standards | ANSI/AAMI ES 60601-1, CAN/CAS C22.2 No.60601.1, IEC/EN60601-1, FDA510(k), FCC Class B, CE, VCCI Class B, KC, ICES-003-B, C-Tick |



WIDE Corporation (Gomae-Dong) 12 Wongomae-Ro, Giheung-Gu, Yongin-Si, Gyeonggi-Do 446-901, Korea Tel: +82-31-218-1600 Fax: +82-31-274-7400 info@widecorp.com **WIDE USA Corporation** 2210 E. Winston Road, Anaheim, CA 92806, USA Tel: +1-714-300-0540 infousa@widecorp.com **WIDE EUROPE B.V.** Panamalaan 3L, 1019 AS Amsterdam, The Netherlands Tel: +31-20-311-9797 Fax: +31-20-311-9790 infoeu@widecorp.com **WIDE JAPAN** Shinjuku Suzuki Bldg A 1-6-8 Shinjuku, Shinjuku-Ku, Tokyo160-0022, Japan Tel : +81-3-6457-8371 Fax : +81-3-6457-8372 infojapan@widecorp.com **WIDE ASIA** 1350/88-89, Thairong tower, Pattanakarn road, Suanluang, Bangkok 10250, Thailand Tel:+66 2 717 1843 Fax: +66 2 717 1844 infoasia@widecorp.com **Foreseeson Custom Displays, Inc.** 2210 E. Winston Road Anaheim, CA 92806 USA Tel: +1 714 300 0540 Fax: +1-714-300 0546 **Foreseeson GmbH** Industriestrasse 38a 63150 Heusenstamm, Germany Tel: +49 6104 64398 0 Fax: +49 6104 64398 11 sales@fsnmed.eu **Foreseeson UK Ltd.** Unit 2 Kingsmill Business Park Chapel Mill Road Kingston upon Thames Surrey KT1 3GZ Tel: +44 208 546 1047 Fax: +44 208 546 3931 sales@fsnmed.eu | Specifications and features are subject to change without notice. Images shown are for illustrative purpose only. All products names are trademarks or registered trademarks of WIDE Corp. This document is copyrighted. Neither this document, nor any part of it, may be reproduced or copied without written permission of WIDE Corp. | Created in Korea. 2016. 02. Ver 1.0.0@MARKSMARKS.NL10-2016

Pattanakarn road, Suanluang, Bangkok 10250, Thailand Tel:+66 2 717 1843 Fax: +66 2 717 1844 infoasia@widecorp.com **Foreseeson Custom Displays, Inc.** 2210 E. Winston Road Anaheim, CA 92806 USA Tel: +1 714 300 0540 Fax: +1-714-300 0546 **Foreseeson GmbH** Industriestrasse 38a 63150 Heusenstamm, Germany Tel: +49 6104 64398 0 Fax: +49 6104 64398 11 sales@fsnmed.eu **Foreseeson UK Ltd.** Unit 2 Kingsmill Business Park Chapel Mill Road Kingston upon Thames Surrey KT1 3GZ Tel: +44 208 546 1047 Fax: +44 208 546 3931 sales@fsnmed.eu | Specifications and features are subject to change without notice. Images shown are for illustrative purpose only. All products names are trademarks or registered trademarks of WIDE Corp. This document is copyrighted. Neither this document, nor any part of it, may be reproduced or copied without written permission of WIDE Corp. | Created in Korea. 2016. 02. Ver 1.0.0@MARKSMARKS.NL10-2016